

Product Datasheet

Clatherin Light Chain Antibody (LCA/LCB) (orb386040)

Description

Recognizes proteins of 31-44kDa, which are


Species/Host

Mouse

Reactivity


Human

Conjugation

Unconjugated

Tested

IHC-P

Applications


Immunohistochemical
staining of Human
pa...

Immunogen

A recombinant N-terminal protein fragment from human Clathrin light chain was used as the immunogen for the Clathrin Light Chain antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the Clathrin Light Chain antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Optimal dilution of the Clathrin Light Chain antibody should be determined by the researcher.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgG2b, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID

[P09496](#), [P09497](#)

Hazard
Information

This Clathrin Light Chain antibody is available for research use only.

Dilution Range

Immunohistochemistry (FFPE): 0.5-1ug/ml for 30 min at RT

Expiration Date

12 months from date of receipt.