

Product Datasheet

Connexin 32 Antibody / GJB1 (orb386045)

Description

This Ab recognizes a protein of 27-32kDa,


Species/Host

Mouse

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IHC-P

Applications

Immunohistochemical staining of Human li...

Immunogen


A recombinant human protein was used as the immunogen for the Connexin 32 antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the Connexin 32 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).



Immunohistochemical staining of Human to...

Note

For research use only

Application notes

Optimal dilution of the Connexin 32 antibody should be determined by the researcher.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgG

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID

P08034

Hazard
Information

This Connexin 32 antibody is available for research use only.

Dilution Range

Immunohistochemistry (FFPE): 0.1-0.2ug/ml for 30 min at RT

Expiration Date

12 months from date of receipt.