

Product Datasheet

KRT15 Antibody / Cytokeratin 15 (orb386052)

Description

Keratins are a family of intermediate filament

forming...

Species/Host

Mouse

Reactivity

Human, Rat

Conjugation

Unconjugated

Tested

IHC-P

Applications

Immunogen

Recombinant human KRT15 protein was used as the immunogen for the Cytokeratin 15 antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the Cytokeratin 15 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Optimal dilution of the Cytokeratin 15 antibody should be determined by the researcher.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgG2a, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID

P19012

Hazard

Information

This Cytokeratin 15 antibody is available for research use only.

Dilution Range

Immunohistochemistry (FFPE): 0.1-0.2ug/ml for 30 min at RT

Expiration Date

12 months from date of receipt.



Immunohistochemical staining of Human basaloid carcinoma tissue.