

Product Datasheet

Glycoprotein 2 Antibody / GP2 / ZAP75 (orb386071)

Description

GP2 (Glycoprotein 2), also known as ZAP75, is a



m...

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC-P, WB

Applications

Immunohistochemical
staining of Human
pa...

Immunogen

A partial human recombinant protein corresponding to amino acids 35-179 was used as the immunogen for the GP2 antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the GP2 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Optimal dilution of the GP2 antibody should be determined by the researcher.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgG2b, lambda

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID

P55259

Hazard
Information

This GP2 antibody is available for research use only.

Dilution Range

ELISA: 2-4ug/ml (order BSA/azide-free format), Western blot: 1-2ug/ml, Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT

Expiration Date

12 months from date of receipt.