

Product Datasheet

hCG Receptor / LH Receptor Antibody (orb386094)

Description

Recognizes a protein of 78-85 kDa, which is


Species/Host

Mouse

Reactivity


Human

Conjugation

Unconjugated

Tested

ELISA, FACS, WB

Applications


Immunohistochemical
staining of Human
ov...

Immunogen

A partial recombinant protein from the N-terminal region of human Luteinizing hormone / Choriogonadotropin receptor was used as the immunogen for the LHR antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the LHR antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Optimal dilution of the LHR antibody should be determined by the researcher.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgG1, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID

P22888

Hazard
Information

This LHR antibody is available for research use only.

Dilution Range

ELISA (order BSA/sodium azide-free format for coating), Flow cytometry: 1-2ug/million cells, Western blot: 0.5-1ug/ml

Expiration Date

12 months from date of receipt.