

## Product Datasheet

### p63 Antibody / Tumor protein 63 (orb386105)

**Description**

p63 is a homolog of the tumor suppressor p53. It


**Species/Host**

Rabbit

**Reactivity**

Human

**Conjugation**

Unconjugated

**Tested**

IHC-P

**Applications**

Immunohistochemical  
staining of Human  
pr...

**Immunogen**

Recombinant partial human p63 protein was used as the immunogen for this p63 antibody.

**Preservatives**

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

**Storage**

Store the p63 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

**Note**

For research use only

**Application notes**

The optimal dilution of the p63 antibody for each application should be determined by the researcher.

**Formula**

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

**Isotype**

Rabbit IgG

**Clonality**

Polyclonal

**Purity**

Protein A affinity chromatography

**Uniprot ID**

**Q9H3D4**

**Hazard**
**Information**

This p63 antibody is available for research use only.

**Dilution Range**

Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT

**Expiration Date**

12 months from date of receipt.