

## Product Datasheet

### Parathyroid Hormone Antibody / N Terminal (orb386112)

## Description

PTH/Parathyroid hormone is a hormone

regulates the concentrat...

## Species/Host

Rabbit

## Reactivity

Human

## Conjugation

Unconjugated

## Tested

## Applications

IHC-P

## Immunogen

A synthetic peptide from the N-terminal region of human Parathyroid hormone was used as the immunogen for the PTH antibody.

## Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

## Storage

Store the PTH antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

## Note

For research use only

## Application notes

Optimal dilution of the PTH antibody should be determined by the researcher.

## Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

## Isotype

Rabbit IgG

## Clonality

Polyclonal

## Purity

Protein A affinity chromatography

## Uniprot ID

**P01270**

## Hazard

## Information


This PTH antibody is available for research use only.

## Dilution Range

Immunohistochemistry (FFPE): 0.5-1ug/ml for 30 min at RT

## Expiration Date

12 months from date of receipt.



Immunohistochemical staining of Human pa...