

Product Datasheet

Growth Hormone Antibody / GH1 (orb386117)

Description

Growth hormone (GH), also known as



[hGH or HGH] in its huma...

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

Applications

IHC-P

Immunohistochemical
staining of Human
pi...

Immunogen

Amino acids 58-187 from the human protein were used as the immunogen for this Growth Hormone antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the Growth Hormone antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Titration of the Growth Hormone antibody may be needed for optimal performance.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgG2b, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID

P01241

Hazard

Information

This Growth Hormone antibody is available for research use only.

Entrez

2688

Dilution Range

Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT

Expiration Date

12 months from date of receipt.