

Product Datasheet

PDCD1 Antibody / PD-1 / PD1 (orb386120)

Description

PDCD1 (Programmed cell death-1 protein), also



recept...

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P

Applications

Immunohistochemical
staining of Human
to...

Immunogen

A full-length human protein was used as the immunogen for the PD1 antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the PD1 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Titering of the PD1 antibody may be required for optimal performance.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Rabbit IgG

Clonality

Polyclonal

Purity

Protein A affinity chromatography

Uniprot ID

[Q15116](#)

Hazard
Information

This PD1 antibody is available for research use only.

Dilution Range

Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT

Expiration Date

12 months from date of receipt.