

Product Datasheet

N-Cadherin Antibody / CDH2 (orb386127)

Description

Recognizes a protein of ~140kDa, identified as N-

Species/Host

Mouse

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested
Applications

FACS, IF, IHC-P, WB

Immunogen

A human partial recombinant protein corresponding to the N-Cadherin extracellular domain was used as the immunogen for this CDH2 antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the CDH2 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

The concentration stated for each application is a general starting point. Variations in protocols, secondaries and substrates may require the CDH2 antibody to be titrated up or down for optimal performance.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgG1, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID

P19022

Hazard
Information


This CDH2 antibody is available for research use only.

Dilution Range

Flow cytometry: 1-2ug/million cells, Immunofluorescence: 1-2ug/ml, Western blot: 1-2ug/ml, Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT

Expiration Date

12 months from date of receipt.



Western blot analysis of Human brain lys...