

Product Datasheet

Thrombomodulin Antibody / CD141 (orb386142)

Description

It recognizes a protein of 75kDa, identified as


Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested


ELISA, IHC-P, WB

Applications

Immunohistochemical
staining of Human
ce...

Immunogen

Human recombinant partial protein corresponding to amino acids 69-194 was used as the immunogen for this Thrombomodulin antibody.


Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the Thrombomodulin antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).


Immunohistochemical
staining of Human
bl...

Note

For research use only

Application notes

The concentration stated for each application is a general starting point. Variations in protocols, secondaries and substrates may require the Thrombomodulin antibody to be titrated up or down for optimal performance.


Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgG2b, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID

P07204

Hazard
Information

This Thrombomodulin antibody is available for research use only.

Dilution Range

ELISA: 2-4ug/ml (order BSA/azide-free format), Western blot: 1-2ug/ml, Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT

Expiration Date

12 months from date of receipt.