

Product Datasheet

Pan Cytokeratin Antibody Cocktail (Acidic + Basic) (orb386219)

Description

Twenty human keratins are resolved with two-



6.0) sub...

Species/Host

Rabbit

Reactivity

Human


Conjugation

Unconjugated

Tested

Applications

FACS, IF, IHC-P, WB



Immunohistochemical
staining of Human
sk...

Immunogen

Recombinant full-length human KRT76 and KRT77 proteins were used as the immunogen for this pan Cytokeratin antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the pan Cytokeratin antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

The concentration stated for each application is a general starting point. Variations in protocols, secondaries and substrates may require the pan Cytokeratin antibody to be titrated up or down for optimal performance.1. The prediluted format is supplied in a dropper bottle and is optimized for use in IHC. After epitope retrieval step (if required), drip mAb solution onto the tissue section and incubate at RT for 30 min.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Rabbit IgG

Clonality

Polyclonal

Purity

Protein A affinity chromatography

Hazard

Information

This pan Cytokeratin antibody is available for research use only.

Entrez

[3866](#), [1](#), [3852](#), [8](#), [3858](#), [4](#), [5](#), [3](#), [14](#), [3868](#), [3853](#), [3861](#), [16](#), [3856](#), [19](#), [3848](#), [3880](#), [15](#), [3851](#), [6](#), [10](#), [3850](#)

Dilution Range

Flow cytometry: 0.5-1ug/million cells, Immunofluorescence: 1-2ug/ml, Western blot: 0.5-1ug/ml for 2 hours at