

Product Datasheet

CD20 Antibody (orb386222)

Description

CD20 is a non-Ig differentiation antigen of B-cells and



Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FACS, IHC-P

Applications

Immunohistochemical
staining of human
to...

Immunogen

Full length human CD20 protein was used as the immunogen for this recombinant CD20 antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Recombinant CD20 antibody with azide can be stored at 2-8°C. The azide-free format should be aliquoted and stored at -20°C or colder.

Note

For research use only

Application notes

Titering of the recombinant CD20 antibody may be required for optimal performance.1. The prediluted format is supplied in a dropper bottle and is optimized for use in IHC. After epitope retrieval step (if required), drip mAb solution onto the tissue section and incubate at RT for 30 min.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Rabbit IgG, kappa

Clonality

Recombinant

Purity

Protein A affinity chromatography

Uniprot ID

P11836

Hazard

Information

This recombinant CD20 antibody is available for research use only.

Dilution Range

Flow cytometry: 1-2ug/million cells, Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT

Expiration Date

12 months from date of receipt.