

## Product Datasheet

### Rabbit Isotype Control Monoclonal Antibody (orb386238)

## Description

Isotype control antibodies are negative

binding or background ...

## Species/Host

Rabbit

## Conjugation

Unconjugated

## Tested

ELISA, FACS, IF, IHC, IP, WB

## Applications

## Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

## Storage

Store the recombinant rabbit isotype control monoclonal antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

## Note

For research use only

## Application notes

This rabbit isotype control monoclonal antibody shows no reactivity to tested lysates and tissues. Generally, the dilution used is matched to the dilution of the antibody being used to test for the protein of interest.

## Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

## Isotype

Rabbit IgG, kappa

## Clonality

Recombinant

## Purity

Protein A affinity chromatography

## Hazard

## Information

This rabbit isotype control monoclonal antibody is available for research use only.

## Dilution Range

ELISA: assay dependent, Flow cytometry: assay dependent, Immunofluorescence: assay dependent, Immunoprecipitation: assay dependent, Immunohistochemistry (FFPE or Frozen): assay dependent, Western blot: assay dependent

## Expiration Date

12 months from date of receipt.