


Product Datasheet

Ku70 + Ku80 antibody (PE) (orb386264)

Description	Mouse monoclonal antibody against Ku70 + Ku80 conjugated to PE
Species/Host	Mouse
Reactivity	Human
Conjugation	PE
Tested Applications	FACS
Immunogen	A human cell nuclear fraction was used as the immunogen for this Ku70 + Ku80 antibody.
Preservatives	500 ul at 0.1 mg/ml with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide
Storage	Store the Ku70 + Ku80 antibody at 2-8°C. Conjugate is light sensitive, store in the dark.
Note	For research use only
Application notes	Titration of the Ku70 + Ku80 antibody may be required for optimal performance.
Formula	500 ul at 0.1 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG1, kappa
Clonality	Monoclonal
Purity	Protein G purified
Hazard Information	This Ku70 + Ku80 antibody is available for research use only.
Entrez	2547
Dilution Range	Flow cytometry: 5ul/test/million cells in 100ul or 5ul/100ul whole blood
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
K562 cells
U...