

Product Datasheet

SFRP5 antibody (orb386269)

Description	Rabbit polyclonal antibody to SFRP5
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human SFRP5
Preservatives	0.01% NaN3
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	36
Uniprot ID	Q9WU66
Expiration Date	12 months from date of receipt.