

Product Datasheet

ZBTB8B antibody (orb386306)

Description	Rabbit polyclonal antibody to ZBTB8B
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human ZBTB8B
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	54
Uniprot ID	Q8NAP8
Expiration Date	12 months from date of receipt.