

Product Datasheet

MMP2 antibody (orb386349)

Description	Rabbit polyclonal antibody to MMP2
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human MMP2
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	74
Uniprot ID	P08253
Expiration Date	12 months from date of receipt.