

Product Datasheet

FKBP11 antibody (orb386358)

Description	Rabbit polyclonal antibody to FKBP11
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of mouse FKBP11
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	22
Uniprot ID	Q9D1M7
Expiration Date	12 months from date of receipt.