

Product Datasheet

GADD45A antibody (orb386359)

Description	Rabbit polyclonal antibody to GADD45A
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	middle region of human GADD45A
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	18
Uniprot ID	P24522
Expiration Date	12 months from date of receipt.