

Product Datasheet

SSU72 antibody (orb386361)

Description	Rabbit polyclonal antibody to SSU72
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	middle region of human SSU72
Preservatives	0.01% NaN3
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	23
Uniprot ID	Q9NP77
Expiration Date	12 months from date of receipt.