

Product Datasheet

RNF4 antibody (orb386399)

Description	Rabbit polyclonal antibody to RNF4
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	middle region of mouse RNF4
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	22
Uniprot ID	Q9QZS2
Expiration Date	12 months from date of receipt.