

Product Datasheet

MKI67IP antibody (orb386400)

Description	Rabbit polyclonal antibody to MKI67IP
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human MKI67IP
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	34
Uniprot ID	Q9BYG3
Expiration Date	12 months from date of receipt.