

Product Datasheet

RALBP1 antibody (orb386410)

Description	Rabbit polyclonal antibody to RALBP1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human RALBP1
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	76
Uniprot ID	Q15311
Expiration Date	12 months from date of receipt.