

Product Datasheet

CSNK1E antibody (orb386415)

Description	Rabbit polyclonal antibody to CSNK1E
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal regions of human CSNK1E
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	47
Uniprot ID	P49674
Expiration Date	12 months from date of receipt.