

Product Datasheet

POU5F1 antibody (orb386438)

Description	Rabbit polyclonal antibody to POU5F1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human POU5F1
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	39
Uniprot ID	Q01860
Expiration Date	12 months from date of receipt.