

Product Datasheet

E2F8 antibody (orb386445)

Description	Rabbit polyclonal antibody to E2F8
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	N-terminal region of human E2F8
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	94
Uniprot ID	A0AVK6
Expiration Date	12 months from date of receipt.