




Product Datasheet

MERTK antibody (orb38645)

Description	Purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	Recombinant Protein
Target	MERTK
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	110249
Uniprot ID	Q12866
NCBI	NP_006334.2
Dilution Range	WB: 1:1000, IF/ICC: 1:10-50, IHC-P: 1:50-10, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Western blot analysis of MERTK Antibody ...



Formalin-fixed and paraffin-embedded hum...