

Product Datasheet

NFKBIB antibody (orb386464)

Description	Rabbit polyclonal antibody to NFKBIB
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	middle region of human NFKBIB
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	38
Uniprot ID	Q15653
Expiration Date	12 months from date of receipt.