

Product Datasheet

FSCN1 antibody (orb386470)

Description	Rabbit polyclonal antibody to FSCN1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human FSCN1
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	55
Uniprot ID	Q16658
Expiration Date	12 months from date of receipt.