

## Product Datasheet

### ZDHHC4 antibody (orb386485)

|                            |                                      |
|----------------------------|--------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to ZDHHC4 |
| <b>Species/Host</b>        | Rabbit                               |
| <b>Reactivity</b>          | Human, Mouse                         |
| <b>Conjugation</b>         | Unconjugated                         |
| <b>Tested Applications</b> | WB                                   |
| <b>Immunogen</b>           | C-terminal region of human ZDHHC4    |
| <b>Preservatives</b>       | 0.01% NaN <sub>3</sub>               |
| <b>Storage</b>             | -20°C                                |
| <b>Note</b>                | For research use only                |
| <b>Clonality</b>           | Polyclonal                           |
| <b>Purity</b>              | Affinity purification                |
| <b>MW</b>                  | 40                                   |
| <b>Uniprot ID</b>          | <b>Q9NPG8</b>                        |
| <b>Expiration Date</b>     | 12 months from date of receipt.      |