

Product Datasheet

E2F7 antibody (orb386487)

Description	Rabbit polyclonal antibody to E2F7
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	N-terminal region of human E2F7
Preservatives	0.01% NaN3
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	100
Uniprot ID	Q96AV8
Expiration Date	12 months from date of receipt.