

Product Datasheet

BAZ2B antibody (orb386512)

Description	Rabbit polyclonal antibody to BAZ2B
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal regions of human BAZ2B
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	240
Uniprot ID	Q9UIF8
Expiration Date	12 months from date of receipt.