

Product Datasheet

RNF1 antibody (orb386521)

Description	Rabbit polyclonal antibody to RNF1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	middle region of human RNF1
Preservatives	0.01% NaN3
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	42
Uniprot ID	Q06587
Expiration Date	12 months from date of receipt.