

Product Datasheet

ANXA4 antibody (orb386529)

Description	Rabbit polyclonal antibody to ANXA4
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human ANX4
Preservatives	0.01% NaN3
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	36
Uniprot ID	P09525
Expiration Date	12 months from date of receipt.