

Product Datasheet

UNC5A antibody (orb386539)

Description	Rabbit polyclonal antibody to UNC5A
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of mouse UNC5A
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	99
Uniprot ID	Q8K1S4
Expiration Date	12 months from date of receipt.