

Product Datasheet

OSBPL11 antibody (orb386546)

Description	Rabbit polyclonal antibody to OSBPL11
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human OSBPL11
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	84
Uniprot ID	Q9BXB4
Expiration Date	12 months from date of receipt.