


Product Datasheet

DOCK5 antibody (orb386549)

Description	Rabbit polyclonal antibody to DOCK5
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human DOCK5
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	215
Uniprot ID	Q9H7D0
Dilution Range	WB: 1:1,000-1:3,000
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell e...