


Product Datasheet

LIN52 antibody (orb386553)

Description	Rabbit polyclonal antibody to LIN52
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human LIN52
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	13
Uniprot ID	Q52LA3
Dilution Range	WB: 1:1,000-1:3,000
Expiration Date	12 months from date of receipt.



Western blot of total cell extracts from...