

Product Datasheet

MAP3K4 antibody (orb386571)

Description	Rabbit polyclonal antibody to MAP3K4
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human MAP3K4
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	182
Uniprot ID	Q9Y6R4
Expiration Date	12 months from date of receipt.