

Product Datasheet

MAPK8IP1 antibody (orb386582)

Description	Rabbit polyclonal antibody to MAPK8IP1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human MAPK8IP1
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	78
Uniprot ID	Q9UQF2
Expiration Date	12 months from date of receipt.