

Product Datasheet

LRRC7 antibody (orb386590)

Description	Rabbit polyclonal antibody to LRRC7
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	middle region of human LRRC7
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	172
Uniprot ID	Q96NW7
Expiration Date	12 months from date of receipt.