

Product Datasheet

GLI1 antibody (orb386595)

Description	Rabbit polyclonal antibody to GLI1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human GLI1
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	118
Uniprot ID	P08151
Expiration Date	12 months from date of receipt.