

Product Datasheet

SNAI2 antibody (orb386598)

Description	Rabbit polyclonal antibody to SNAI2
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human SNAI2
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	30
Uniprot ID	043623
Expiration Date	12 months from date of receipt.