

Product Datasheet

SMC6 antibody (orb386618)

Description	Rabbit polyclonal antibody to SMC6
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human SMC6
Preservatives	0.01% NaN3
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	126
Uniprot ID	Q96SB8
Expiration Date	12 months from date of receipt.