

Product Datasheet

ZBTB1 antibody (orb386619)

Description	Rabbit polyclonal antibody to ZBTB1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human ZBTB1
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	82
Uniprot ID	Q9Y2K1
Expiration Date	12 months from date of receipt.