

Product Datasheet

ZDHHC1 antibody (orb386631)

Description	Rabbit polyclonal antibody to ZDHHC1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human ZDHHC1
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	55
Uniprot ID	Q8WTX9
Expiration Date	12 months from date of receipt.