

Product Datasheet

DESMIN antibody (orb386643)

Description	Rabbit polyclonal antibody to DESMIN
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human DESMIN
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	53.5
Uniprot ID	P17661
Expiration Date	12 months from date of receipt.