

Product Datasheet

SIRT4 antibody (orb386650)

Description	Rabbit polyclonal antibody to SIRT4
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human SIRT4
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	35
Uniprot ID	Q9Y6E7
Expiration Date	12 months from date of receipt.