

Product Datasheet

ZSCAN16 antibody (orb386651)

Description	Rabbit polyclonal antibody to ZSCAN16
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	middle region of human ZSCAN16
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	41
Uniprot ID	Q9H4T2
Expiration Date	12 months from date of receipt.